

Flight Case – Portable X-Ray Kit

Ready for any situation – this flight case keeps all your dental X-ray imaging equipment together and secure.

The portable Flight Case from INLINE Systems enables you to carry all the required dental imaging and diagnostic equipment for any examination in a single transport case. The case is fitted with a telescopic handle and castors for ease of use, as well as with a custom foam insert to guarantee safe transport.

Whether the equipment is being transported within a clinic or hospital, or across the country in the back of a 4WD, the Flight Case will ensure that everything is in perfect condition and charged when you need it.

**Portable
imaging
doesn't get
any easier.**

**exclusive
distributor**



Flight Case – Portable X-Ray Kit



Keep your equipment charged and ready for use.

Take quality images with a truly portable hand-held x-ray.

All the equipment is battery powered and can also be used with mains power. When inserted into the case, each piece of equipment is connected to a central charging circuit with a single power connect for ease of use. Simply plug the power lead into the mains socket and each piece of equipment in the case will begin charging.

Making the equipment in the Flight Case truly portable, the configuration includes the battery powered NOMAD Pro 2™ handheld X-ray unit. This system lets you gain all the reliability and efficiency advantages that every busy practitioner needs and appreciates. It provides a powerful freedom in workflow, improving dental radiography speed, convenience, and efficiency thanks to its lighter weight and cordless design.

X-ray images are transmitted back to the toughened DuraBook PC and appear on the screen in a matter of seconds after the exposure.

Features & benefits:

- Portable
- Durable
- Compact
- Safe transport

Kit Content

- Transport case with custom foam insert
- NOMAD Pro 2™ handheld x-ray unit
- Two NiCad batteries
- Mains charger unit for batteries
- Digital x-ray system
- Digital x-ray sensor (size 1 or 2 or both)
- USB2 intraoral camera
- DuraBook toughened notebook PC, Win-XP
- DentalEye imaging software
- Canon photo printer, battery operated
- Document satchel in lid of case

KAVO

Dental Excellence

ISO 13485 Quality Certified

For further information and a demonstration, contact **Inline** at info@inline.com.au

1300 033 733

inline.com.au

 **Inline**
IMAGING TECHNOLOGY